

FRANCK MULLER
GENEVE

FRANCK MULLER
GENEVE

Vanguard Rose
Skeleton



V 32 S6 SQT ROSE (RG)

5N



Rose Collection

The Vanguard™ Rose Skeleton combines beautiful design with ingenious development. The technicality of the collection lies in the research aimed at making the roses integral to the mechanism with an openwork dial. Here, the flowers are not mere ornaments placed on the movement; instead, they are cut from the bridges and the bottom plate, connecting the mechanical movements together. They serve as an aesthetic decoration essential for this beautiful timepiece.

Reference V 32 S6 SQT ROSE (RG)

Functions Hours, Minutes, Seconds at 6 o'clock

Movement MVT FM 1540-VS11
 In-House Movement
 Manually wound movement
 4 days power reserve
 167 components
 Width: 26,75 mm. Length: 32 mm. Thickness: 5,9 mm
 Balance wheel frequency set to 21,600 vibrations per hour
 Skeleton movement
 Two barrels

Movement Decoration Mirror polishing on the regulator assembly
 Beveling on the mainplate, bridges and regulator assembly
 Flank drawing on the bridges and mainplate
 Diamond polishing on the sinks
 Circular satin brushing on the wheels, drums and barrel covers
 Satin brushing on the mainplate and bridges
 Enamel painting on the hearts

Instructions Dual-Position Crown: 1. Winding. 2. Setting the time

Case Vanguard case
 Water resistant up to 30 meters
 Width: 32 mm. Length: 42.3 mm. Thickness: 9.9 mm
 18k rose gold case
 Hand polished
 Red inserts on both sides of the case
 Sapphire crystal

Bracelet Alligator strap

Buckle 18k rose gold buckle

